

Computer minimal requirements for PRIMO system

Minimal configuration :

- Workstation DELL PRECISION type T3620
- OS : Windows 7 and 8 Professionnal 64-bits
- Processeur : Intel Xéon Processeur E3-1226 v3 (Quad core, 3,30 GHz Turbo, 8 Mo) or equivalent
- RAM : 16 Go 1600 MHz DDR3 ad minima
- Ad minima 1 Port Display or 1 port HDMI available for the DMD
- 2 available PCI slots may be needed:
 - 1 for the USB3 card of Hamamatsu Flash4LT camera
 - 1 for the Marzhauser stage PCI card
- Monitor full HD with 1920 x 1080 resolution
- Hard disk : 500 Go SATA 7200 tr/min 3,5"
- DVD player: 16x DVD+/-RW Player